GAP ENDO-EXO MEDULLARY SYSTEM

ENGLISH

Device Description

The GAP Endo-Exo Medullary System is used for the treatment of fractures or correction of deformities in the femur, tibia and humerus of pediatric patients (child and adolescent). This fixation device consists of an intramedullary nail linked to a plate via a lag and/or mechanical screws creating a combined Endomedullary/Exomedullary osteosynthesis device. The novel approach of osteosynthesis intends to create a load sharing system between the nail and plate with the objective of limiting the risk of stress fractures and improving the implant stability. The plates have two main functions: a) to give lateral support to weak lateral cortex of an osteoporotic bone avoiding concentrated stresses at the screw/head/bone interface; b) to lock the lag screws when treating femoral neck fractures or to lock the k-wires when treating coxa vara or coxa valga using a sub lesser trochanter osteotomy. The plates are part of the system and not intended to be used as stand-alone fracture plates.

Washers can be used under the cortical screw heads to increase the contact surface between the screw and the host bone, reducing the possibility of stress fracture in the bone. The intramedullary nail's fixation to the bone is assured by a conical thread at its head and by one or several transverse cortical screws at its tip.

Intended Use

The GAP Endo-Exo Medullary System is indicated as a temporary implant to ensure alignment, stabilization and fixation of: pathological long bones that have been surgically prepared (osteotomy) for correction of deformities or fractures caused by trauma or disease.

The GAP System is used for pediatric patients (child and adolescent) with skeletal dysplasias.

It can be used to correct the following conditions:

- Diaphyseal fracture of the femur, tibia and humerus
- Fractures of the femoral neck
- Subtrochanteric, intertrochanteric and combination fractures
- Correction of deformities (OI, Coxa vara, Coxa valga)
- Nonunions and malunions

Contraindications

- Active or suspected latent infection or marked local inflammation in or about the affected area.
- Osteoporosis, insufficient quality or quantity of bone/soft tissue
- Compromised vascularity inhibiting adequate blood supply to the operative site.
- Documented or suspected material sensitivity.
- Patients with abnormal neurological or mental conditions
- Sepsis
- Other medical or surgical conditions which would preclude the potential benefit of surgery.

Surgeons should warn patients about these contraindications and limitations when appropriate.

Adverse Effects

- Pain, discomfort or abnormal sensations due to the presence of the device.
- Limb shortening or residual deformity with nonunion or malunion
- Metal sensitivity and/or allergic reaction to a foreign body.
- Nerve damage due to the surgical trauma
- Bone resorption due to stress shielding.
- Postoperative bone fracture and pain.
- Infection, both deep and superficial
- Unrecognized joint penetration
- Inadequate healing
- Necrosis of bone

Warnings

- Pega Medical advises against the use of another manufacturer's component with any Pega Medical component. Any such use will negate the responsibility of Pega Medical for the performance of the resulting mix.
- Implants are single use items. Please note that single use devices (SUD) that come into contact with human blood or tissue should not be reused and should be returned to the manufacturer or disposed of properly.
- Metal implants should never be reused. Although appearing undamaged, the device may have small defects or internal stresses that may cause implant failure.
- Correct implant handling is extremely important. Avoid contouring of metallic implants. Discard all damaged or mishandled implants, or return them to the manufacturer for proper disposal.
- Continuous screening with an image intensifier (fluoroscopy) during guide wire insertion and whenever cannulated instruments are advanced over a guide wire is recommended to prevent unintended guide wire advancement and penetration into the surrounding tissues.
- Improper insertion of the device during implantation can increase the possibility of loosening or migration.
- Selecting the largest diameter implant that is appropriate for the medullary canal of the host bone as well as proper positioning and insertion of the implant are crucial to mitigate the risk of implant failure.
- Device breakage or damage can occur when implant is subjected to increased loading associated with delayed union, non-union, or incomplete healing. Proper consolidation should be observed prior to full weight bearing.
- Plate Bending: The plate should not be excessively or repeatedly bent. The plate should not be reverse bent in the same location. Use care to ensure that plate is not scratched or notched during the bending process.
- Contouring and bending of an implant may reduce its fatigue strength causing failure under load.
- Screws and plates included in the GAP Endo-Exo Medullary system can only be used with the GAP Nail. The plates included in the GAP system are not standalone osteosynthesis plates.
- A minimum of two Cortical Screws must be used for distal fixation of the Nail.
- Implant System can only be used for patients weighing 60 kg and under or as indicated in the table below.
- For fractures or osteotomies below the lesser trochanter combined with Lag Screw use, the following Lag Screw and weight limitations should be observed:

| Nail Size (Ø) | Max. Allowable Lag Screw Length | Max. Patient Weight |
|---------------|---------------------------------|---------------------|
| 4.8 | 50 mm | 40 Kg |
| 5.6 | 70 mm | 40 Kg |
| 6.4 | 80 mm | 50 Kg |
| 7.2 and above | No Limit | 60 Kg |

MRl Safety Information

The GAP Endo-Exo Medullary System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the GAP Endo-Exo Medullary System in the MR environment is unknown.

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

Surgical Technique

Standard surgical technique for the placement of the intramedullary nails is recommended. Pre-operative procedures, knowledge of applicable surgical techniques, proper patient selection and correct implant placement are all equally important for the successful use of this product. The surgical technique manual details every step and should be carefully followed.

Device Lifetime and Retrieval

Removal or replacement of the implant is recommended subsequent to normal follow-up after the bone has consolidated and the deformity correction has been achieved. Routine removal of internal fixation devices may reduce the occurrence of symptomatic complications of implant breakage, implant loosening and implant related pain. In addition, if removal is favorable, Pega Medical recommends the retrieval of implants in order to avoid bone reduction and weakening, particularly in young and active patients. Ensure that bone consolidation is complete prior to the removal of the device. Although the final decision to retrieve the implants falls on the surgeon, a maximum Device Lifetime of 5 years for the implant has been defined to ensure material stability. The Surgical Technique manual details retrieval steps and should be carefully followed.

Cleaning and Sterilization Instructions for Implant Components

All implants are provided clean, but are NON-STERILE when shipped from Pega Medical. The instructions below should be followed for sterilizing items supplied non-sterile. Apply a standard cleaning protocol that is approved by the hospital before implant sterilization. All metallic implants can be steam sterilized following the instructions and parameters listed below:

Note: For USA only: Sterilization wraps, pouches, indicators and sterilization trays should be FDA-cleared for the sterilization cycle parameters.

- Implant components of the GAP System should be sterilized using sterilization pouches.
- Devices should be dry before packaged for sterilization

| Method | Steam |
|---------------------|---------------|
| Sterilization type | Prevacuum |
| Minimal temperature | 270°F (132°C) |
| Minimal cycle time | 4 minutes |
| Minimal drying time | 30 minutes |

Warning: Do not stack trays during sterilization

Other sterilization methods and cycles may also be suitable. However, validation of any alternative method using appropriate laboratory techniques is advised.

Cleaning, Sterilization and Re-sterilization Instructions for Instruments

Reusable instruments must be cleaned and sterilized prior to every use. The instrument tray and instruments of the GAP system should be sterilized wrapped in two layers of 1-ply polypropylene wrap using sequential wrapping techniques. Please refer to document entitled: “Guidance for Instrument Care” for further information and instructions regarding cleaning, sterilization and re-sterilization of instruments.

Notice to the User and/or Patient

If you experience or are aware of any serious incidents that have occurred in relation to the device, please report them to Pega Medical Inc. at feedback@pegamedical.com. Additionally, kindly notify the competent authority of the Member State in which the user and/or patient is established. For Switzerland, please report them to materiovigilance@swissmedic.ch.

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