**Device description**

The Hinge Pediatric Plating System is an articulated plate used for guided hard tissue growth and deformity correction. It includes one plate and two screws. By fixing the plate by means of the screws at each side of the the physaeal plate, axial growth can be controlled and angular deformity corrected.

Plates and screws are manufactured in Medical grade Stainless Steel (316L as per ASTM F138).

The plate is provided in three different sizes with distances between the anchoring points of the screws of 12mm, 16mm and 20mm. The 4.5 mm diameter screws are provided in three different lengths of 25mm, 30mm and 35mm.

**Intended use**

This implant system is indicated as a temporary implant to aid in the correction of the angle of growth of long bones by inhibiting longitudinal growth of the physis in pediatric patients. It can be used to correct the following conditions:

- Femur and tibia: varus, valgus, flexion or extension deformities of the knee.
- Humerus: valgus or varus deformities of the elbow.
- Radius and Ulna: flexion or extension deformities of the wrist.
- Ankle: varus, valgus or plantar flexion deformities of the ankle.

**Contraindications**

- Physaeal bar – unresectable
- Skeletal maturity
- Physiologic genu varum/valgum
- Osteoporosis, insufficient quality of quantity of bone/soft tissue
- Material sensitivity
- Sepsis
- Patients with abnormal neurological or mental conditions

**Adverse Effects**

- Overcorrection or under-correction of the angular deformity if the device is not removed at the appropriate time.
- Device breakage or damage due to increased loading associated with the treatment if the device is not removed at the appropriate follow up time.
- Possibility of loosening or migration of the implant due to improper insertion of the device during implantation.
- Pain, discomfort, or abnormal sensations due to the presence of the device.
- Bending, fracture or migration of the implant
- Metal sensibility and/or allergic reaction to a foreign body.
- Nerve damage due to the surgical trauma
- Postoperative bone fracture and pain.
- Infection; both deep and superficial.
- Unrecognized joint penetration.
- Inadequate healing.
- Necrosis of bone.
- Haematoma.

**Warnings**

- Pega Medical advises against the use of another manufacturer’s component(s) with any Pega Medical component(s). Any such use will negate the responsibility of Pega Medical for the performance of the resulting mix.
- Implants are single use items. Please note that single use devices (SUD) that come into contact with human blood or tissue should not be reused and should be returned to the manufacturer or disposed of properly.
- Metal implants should never be reused. Although appearing undamaged, the device may have small defects or internal stresses that may cause implant failure.
- Correct implant handling is extremely important. Avoid contouring of metallic implants. Discard all damaged or mishandled implants, or return them to the manufacturer for proper disposal.
- Continuous screening with an image intensifier (fluoroscopy) during guide wire insertion and whenever cannulated instruments are advanced over a guide wire is recommended to prevent unintended guide wire advancement and penetration into the surrounding tissues.
- Improper insertion of the device during implantation can increase the possibility of loosening or migration.

**MRI Safety Information**

The Hinge Pediatric Plating System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Hinge Pediatric Plating System in the MR environment is unknown.

**CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Surgical technique and follow-up**

Precise placement of the Hinge Pediatric Plating System is achieved by means of device-specific instrumentation that allows positioning of guidance pins on a template. With the template in place, predrilling of the screw holes can be done under image intensifier over guidance pins with cannulated drills. An implant holder completes the precision instrumentation. The surgical technique manual should be carefully followed.

Avoid any damage of the periosteum during insertion of removal of the device. A physaeal bar could be formed and generate further angular deformity.

It is recommended a three months follow-up of patients during the implantation period and an annual follow-up after removal until maturity to assure maintenance of correction.

Slight overcorrection might be necessary, especially in young patients, to compensate rebound deformity after removal.

**Device Lifetime and Retrieval**

This device is used to correct an angular deformation. Once that deformation has been corrected, or slightly over-corrected to account for rebounding after removal of the device, the implant should be removed to prevent adverse over-correction from occurring.

Typically, the Hinge Pediatric Plating System will be implanted for a number of months; occasionally, it will remain implanted for up to three years. It is recommended to follow-up with patients at least every three months during the implantation period, and annually after removal until bone maturity is reached. The final decision to recover the implant falls to the surgeon. It is based on careful follow-up of the patient’s overall condition and the bone’s healing.

**Cleaning, sterilization and re-sterilization instructions for Implants**

All implants are provided clean, but are NON-STERILE when shipped from Pega Medical. The instructions below should be followed for sterilizing items supplied non-sterile. Apply a standard cleaning protocol that is approved by the hospital before implant sterilization. All metallic implants can be steam sterilized following the instructions and parameters listed below:

**Note: For USA only: sterilization wraps, pouches, indicators, and sterilization trays should be FDA-cleared for the sterilization cycle parameters.**

- Implant components of the Hinge Plate System should be sterilized using sterilization pouches
- Devices should be dry before packaging for sterilization

Method	Steam
Sterilization type	Prevacuum
Preconditionning pulse	3
Minimal temperature	270°F (132°C)
Minimal cycle time	4 minutes
Minimal drying time	20 minutes

**Warning: Do not stack trays during sterilization**

Other sterilization methods and cycles may also be suitable. However, validation of any alternative method using appropriate laboratory techniques is advised.

**Cleaning, sterilization and re-sterilization instructions for Instruments**

Reusable instruments must be cleaned and sterilized prior to every use. The instrument tray and instruments of the Hinge Plate System should be sterilized, wrapped in two layers of 1-ply polypropylene wrap, using sequential wrapping techniques.

Please refer to document entitled: “Guidance for Instrument Care” for further information and instructions regarding cleaning, sterilization, and re-sterilization of instruments.

**Notice to the User and/or Patient**

If you experience or are aware of any serious incidents that have occurred in relation to the device, please report them to Pega Medical Inc. at [feedback@pegamedical.com](mailto:feedback@pegamedical.com). Additionally, kindly notify the competent authority of the Member State in which the user and/or patient is established. For Switzerland, please report them to [materiovigilance@swissmedic.ch](mailto:materiovigilance@swissmedic.ch).

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