Device Description

The Locking Pediatric Osteotomy Plate (LolliPOP™) system is a modular hip blade–plate system intended for stable fixation of valgus, varus, derotation, flexion and extension, of proximal femoral osteotomies (PFO) and fractures in the pediatric population. PFOs are widely performed reconstructive surgeries in children with hip deformities, such as Coxa Valga, Coxa Vara and other congenital deformities.

The implants are made of medical grade 316L Stainless Steel (as per ASTM F138) and is offered in four sizes: infant, child, adolescent and HD adolescent. All required instrumentation is provided in the instrument case.

Indications for Use

The Locking Pediatric Osteotomy Plate (LolliPOP™) system is a modular hip blade–plate system intended for stable fixation of valgus, varus, derotation, flexion and extension of proximal femoral osteotomies (PFO) in the pediatric (infant, child and adolescent) population.

Intended uses include the following:

- Inter and subtrochanteric valgus osteotomies
- Inter and subtrochanteric varus osteotomies
- Inter and subtrochanteric derotation osteotomies
- Inter and subtrochanteric flexion and extension osteotomies
- Inter and subtrochanteric fractures

CAUTION: Federal law (USA) and most other countries' laws restrict this device to sale by or on the order of a physician.

NOTE: This product is marketed for the specific indications described in the labeling. The use of this product for other than its intended purpose(s) is either contraindicated (see CONTRAINDICATIONS) or is without evidence to support the safety and effectiveness of such use. For the information of individuals and institutions contemplating the use of this product for other than labeled indication (i.e., off-label use), such use may be experimental and may be the subject of restrictions under applicable laws and regulations.

Contraindications:

- Do not use in any situation that is not described in the INDICATIONS section of this insert.
- Do not use in patients without adequate tissues coverage at the implant site.
- Do not use in patients with known metal sensitivity or intolerance.
- Do not allow physical contact of the plate, blade and screws of this system with any implant made of different or incompatible metals.
- Do not use in anatomical location in which the device would interfere with other critical structures such as nerves, blood vessels, or other vital structures.
- Do not use in presence of documented infection, suspected latent infection, or marked local inflammation in and about the affected area.
- Do not use in presence of compromised vascularity inhibiting adequate blood supply to the operative site.
- Do not use in presence of severe muscular, neural or vascular diseases
- Do not use in the presence of severe osteopenia and/or osteoporosis, insufficient quality or quantity of bone/soft tissue, or in the presence of marked or rapid bone absorption, metabolic bone disease, sepsis, cancer, or any other tumor-like condition of the bone which may compromise fixation.
- Do not use this system in any medical or surgical situation that would preclude the benefit of surgery such as undiagnosed infection, end stage malignant disease, or other unexplained diseases.
- In skeletally immature individuals, long bone epiphysis or trochanteric epiphysis should not be violated by the device. Cessation of growth may take place.

Adverse Effects:

The risks associated with this device are the same as with any metallic internal fixation device. These include, but are not limited to the following:

- Pain, discomfort or abnormal sensations due to the presence of the device.

- Irritation or inflammation of surrounding soft tissue or skin over implant if coverage is insufficient
- Infection, both deep and superficial
- Metal sensibility and/or allergic reaction to a foreign body.
- Nerve damage resulting in temporary or permanent loss of neurological function.
- Possible neurovascular injury
- Possible blood circulation or vessel damage, or avascular necrosis (AVN)
- Cessation of growth of the operated portion of bone
- Bony formation around implant making removal difficult or impossible
- Bone resorption due to stress shielding leading to peri-implant fracture
- Persistent instability
- Nonunion or delayed union of bone fracture or bony fusion
- Limb shortening or residual deformity with nonunion or malunion
- Migration resulting in injury to soft tissue, visceral organs, or unrecognized joint penetration
- Cutting of osteoporotic bone with loss of reduction
- Corrosion of implants
- Additional injuries can occur as a result of postoperative trauma.
- Device breakage, loosening, and/or migration can occur as a result of early weight bearing or muscle activity.
- It may be necessary to perform additional surgery in order to correct adverse effects or reactions which may not be related to the actual system.

These adverse effects include adverse effects that are important considerations for metallic internal fixation devices. These risks and general surgical risks should be explained to the patient prior to surgery.

Warnings

- Implants are single use items. Please note that single use device (SUD) which comes into contact with human blood or tissue should not be re-used and should be returned to the manufacturer or properly disposed.
- Metal implants should never be re-implanted. Although appearing undamaged, the device may have small defects or internal stresses which may eventually lead to implant failure.
- Correct implant handling is extremely important in order to avoid notching or scratching of the surface of the device. It is also recommended that excessive contouring and bending of an implant be avoided. Contouring or bending of the device may reduce its fatigue strength causing failure under load. Discard all damaged or mishandled implants or return to manufacturer for proper disposal.
- The blade is sharp, handle with care in order to avoid harm to the patient or surgical staff.
- Failure to use largest possible components or improper positioning/insertion of the device during implantation can increase the possibility of migration, loosening, bending, cracking, or fracture of the device or bone, or both. However, the following weight limitations should be observed:

Implant Size	Max Patient Weight
Infant	15 Kg
Child	35 Kg
Adolescent & HD adolescent	55 Kg

- Continuous screening with an image intensifier (fluoroscopy) during guide wire insertion and whenever cannulated instruments or the blades are advanced over a guide wire is recommended to prevent unintended guide wire advancement and penetration into the surrounding tissues.
- This system is designed to be adjunct in the fixation of fractures, dislocations, and other bone fusion procedures, as a temporary implant. No implant can withstand body loads indefinitely without the healing of bone and-or ligaments, the system should be used to augment bone fusion and secure fractures until union occurs. Device

- loosening, disassembly, and/or breakage can occur when implant is subjected to increased loading associated with delayed union, non-union, or incomplete healing.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management to avoid re-fracture or recurrent deformity.
- The patient's mobility should be restricted at the region of the osteotomy or fracture to allow bony union. If a nonunion develops, the implants should be removed. If a solid fusion of bone does not occur, the site should be immobilized until solid bony fusion can be achieved. Failure to immobilize a delayed or nonunion of bone will result in excessive and repeated stresses which are transmitted by the body to any temporary internal fixation device prior to healing of the fracture. Due to normal metal fatigue these stresses can cause eventual bending or breakage of the device.
- Screws, blades and plates included in the Locking Pediatric Osteotomy Plate (LolliPOP™) system can only be used within this system. They are not stand-alone osteosynthesis components.
- Pega Medical advises against the use of another manufacturer's component with any Pega Medical component. Any such use will negate the responsibility of Pega Medical for the performance of the resulting mix.
- Care should be taken not to cut through surgical gloves when handling any sharp-edged implants and instruments, and to take into account the risk of infection if a cut appears.

The surgeon should be aware and the patient informed of the following information and limitation.

- Compliance of the patient may affect the results of the fixation.
- Incidence of nonunion is increased in patient known to smoke.
- Patients should be warned to avoid any sudden change in position, strenuous activity, or falls. To achieve a successful union, the patient should not be exposed to mechanical vibrations, whether intrinsic or extrinsic, that may lead to loosening of the device. The patient should be warned of this possibility and instructed to restrict physical activities especially those causing any type of mechanical stress on the area that is being secured by the system. The patient should avoid any type of sport activities or strenuous work during the postoperative healing period.

Complications and/or failure are more likely to occur in:

- Physically active patients
- Debilitated patients or patients unable to follow instruction or use weight supporting devices
- Patients that suddenly change position, fall, or are exposed to mechanical vibrations.

MRI Safety Information:

The Locking Pediatric Osteotomy Plate (LolliPOP™) system has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Locking Pediatric Osteotomy Plate (LolliPOP™) system in the MR environment is unknown.

Surgical Technique

Pre-operative procedures, knowledge of applicable surgical techniques, proper patient selection and correct placement of the implants are all equally important for the successful use of the Locking Pediatric Osteotomy Plate (LolliPOP™) system. The *Surgical Technique* manual details every step and should be carefully followed.

Device Lifetime and Retrieval

Removal of the implant is recommended subsequent to normal follow-up after the bone has consolidated and the deformity correction has been achieved. Routine removal of internal fixation devices may reduce the occurrence of symptomatic complications of implant breakage, implant loosening and implant related pain. In addition, if removal is favorable, Pega Medical recommends the retrieval of implants in order to avoid bone reduction and weakening, particularly in young and active patients. Ensure that consolidation is complete prior to the removal of the device. Although the final decision to recover the implants falls on the surgeon, a maximum Device Lifetime of 5 years for

the implant has been defined to ensure material stability. The Surgical Technique manual details retrieval steps and should be carefully followed.

Cleaning and Sterilization Instructions for Implant Components

All implants are provided clean, but are NON-STERILE when shipped from Pega Medical. The instructions below should be followed for sterilizing items supplied non-sterile. Apply a standard cleaning protocol that is approved by the hospital before implant sterilization. All metallic implants can be steam sterilized following the instructions and parameters listed below:

Note: For USA only: Sterilization wraps, pouches, indicators and sterilization trays should be FDA cleared for the sterilization cycle parameters.

- **Implants** should be sterilized using sterilization pouches
- **Implant trays** should be sterilized wrapped in two layer of 1-ply polypropylene wrap using sequential wrapping techniques.
- **Devices should be dry before being packaged for sterilization**

Method	Steam
Sterilization type	Prevacuum
Minimal temperature	270°F (132°C)
Minimal cycle time	4 minutes
Minimal drying time	60 minutes

Warning: Do not stack trays during sterilization

Other sterilization methods and cycles may also be suitable. However, validation of any alternative method using appropriate laboratory techniques is advised.

Cleaning, Sterilization and Re-sterilization Instructions for Instruments

Reusable instruments must be cleaned and sterilized prior to every use. The instrument tray and instruments of the Locking Pediatric Osteotomy Plate (LolliPOP™) system should be sterilized wrapped in two layers of 1-ply polypropylene wrap using sequential wrapping techniques.

Please refer to document entitled "*Guidance for Instrument Care*" for further information and instructions regarding cleaning, sterilization and re-sterilization of instruments.

For any further questions, please contact Pega Medical

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www.medicalex.info

CE0413

Manufacturer / Fabricante / Fabricant

Caution: law prohibits dispensing without prescription

See instruction for use / Ver instrucciones de uso / Lire les instructions avant usage

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Authorized representative / Representante Autorizado / Représentant Autorisé

Batch number / Número de lote / Numéro de lot

The Locking Pediatric
Osteotomy Plate System™

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© 2021 Pega Medical Inc.